

Technical specifications

Model m.u.	Laser	Wavelength mm	Pulse Width * ms	Max Repetition Rate Hz	Max Fluence J/cm2	Max energy [J]	Max power at 3 Hz [W]
Alex Shine ^{HT}	Alexandrite	755	1-300	2	15@20mm	50	120
Alex Shine Turbo ^{HT}	Alexandrite	755	1-300	2	24@20mm	80	170
Yag Shine ^{HT}	Nd:YAG	1064 long p	1-300	7	32@20mm	100	220
		1064 short p	0,3-1,5	9	20@8mm	10	30
Dual Shine Blended ^{HT}	Alexandrite	755	1-300	3	21@20mm	66	150
	Nd:YAG	1064 long p	1-300	7	32@20mm	100	220
		1064 short p	0,3-1,5	9	20@8mm	10	30
	Alex + Nd:YAG blended		1-300	2	25@20mm	85	180
Cooling system	All Systems						

* automatic slicing of the pulse

Weight & Dimensions	75 kg. 56W x 88L x 98H cm
Electrical requirements	220-230V AC; 16 A; 50-60 Hz

Handpiece

Standard Spot Size	1-5 mm, 8 mm, 14 mm, 18 mm, 20 mm
Optional Spot Size	10 mm, 12 mm, 16 mm, 22 mm, 24 mm, 26 mm



New hummer handpiece up to 26mm with automatic spot detector



Connection between cooling and handpiece



MIA LASER SRL
Via D. Manin 60
30174 Mestre - Venezia
ITALY

P.IVA IT04289660278
mialaser@mialaser.biz
www.mialaser.biz

not ordinary people

The Blue Series HT lasers are the most efficient in the market and the cutting edge of the aesthetic medicine. Based only on high power technology, these lasers are designed to grant the top of reliability, safety and efficacy to every customer.

The fastest and most effective laser in the market

The series includes:

- Alex Shine HT
- Alex Shine Turbo HT
- Yag Shine HT
- Dual Shine Blended HT




Upgradable to:

- Plus (with IPL integrated)
- Air Cooling

Covering the following applications:

- Hirsutism
- Vascular Lesion
- Telangiectasias
- Spider Angiomas
- Spider Veins
- Red Superficial Veins
- Blue Superficial Veins (less than 4mm in diameter)

FEATURES AT 2Hz @ 20mm SPOT SIZE



Photos courtesy of F.Iacoponi, MD

Customer Review

"I'm using MiaLaser's high end Dual Shine Blended HT since 6 months and I never thought to find a such amazing system, delivering top power and fast repetition. Thanks to the 26 mm spot size at 3 HZ it takes only 45 minutes for a full body treatment. But the most impressive thing is the final result, due to the extreme technology of the machine, that allows to achieve an excellent reduction of hair grow in less sessions and the pain is almost zero.

Dr. Momammed Al Bannai , Cavallo Clinic UAE