

MIA FRAXEL ULTRA PULSED CO2 LASER 30-40 W

Available special kit for :
VAGINAL REJUVENATION



Produced by our partner

SOURCE	Laser CO2
SPECTRUM OF EMISSIONS	10.600 nm
EMISSION MODE	Continous, pulsed, single pulsed, pulse frequency, ultrapulse
DIVERGENCY	< 1,5 mrad
POWER	0.5-30/40 W
SUPER PULSE POWER	500 watt
PULSE FREQUENCY	1 to 500 Hz
PULSE DURATION	10-1000 msec (optional 4000)
ACCESSORIES	Scanner, Surgical Handpiece, Vaginal Rejuvenation Kit

APPLICATIONS

- Dermatology
- Gynecology
- Proctology
- Vet

- Vaginal Rejuvenation
- Urology
- Neurosurgery
- Lentigo removal

SCANNER

Microspotsize
Spot Density
Burst Mode

200 μ adjustable

0,2 – 2 mm

ton 10-1000msec
toff 50-1500msec

Frame

square; square with hole;
triangle;rectangle,hezagon

